Loaner / Rental Equipment Checklist

Make / Model:	Serial Number:
Warehouse Checklist (check box if item reviewed is clean and operational)	
☐ Upholstery and Arm Pads (check for cracks and wear) ☐ Armrests (check / adjust)	
☐ Seating System (check brackets)	
Pneumatic Tires (check wear / pressure)	
Solid Tires (check wear / tightness) Wheels (spokes / mags)	
☐ Hand Rims and Mounting Hardware ☐ Wheel Locks (check for tightness / function	
☐ Wheel Locks (check for tightness / function☐ Casters (check for wear)	''
Bearings (check caster / stem condition and	d adjust)
Front Riggings (check foot plates / heel loops)	
Frame (check side frame / cross braces / side posts)	
Folding (ensure chair folds / unfolds properly)	
☐ Tips and Grips (check for wear)	''',
Pulleys, Belts, Clutches (check for wear)	
☐ Cables (check for visual wear)	
☐ Tilt (ensure tilt is tested with weight in it)	
Test Batteries (result:)
☐ Joystick / Controller (check for secure mounting and function)	
□ Charger	,
☐ Reset Switch (present and operational)	
Anti-tippers (present / functional)	
☐ General (ensure chair is fully operational in all directions / modes, and NSM sticker is attached)	
Checked by:	Date:
Supervisor QA:	Date:
Delivery Technician Checklist (check box if item reviewed is clean and operational)	
☐ Clean and Fully Operational (NSM sticker attached)	
☐ Cushion and Back (check for wear and proper function / ensure you have pump if applicable)	
☐ Seat to Floor Height:	
Delivered by:	Datos